

Corning® U-Shaped Cell Culture Flask CORNING

SAME:

- Surface treatment
- cGMP manufacturing
- Cell growth area
- Stable stacking design
- Large, clear viewing surface
- High standards of quality (100% leak tested)
- Nonpyrogenic (<0.1 EU/mL)
- BSE/TSE compliant
- Canted neck for easier pouring and pipet access

Corning 75 cm² Cell Culture U-shaped Flasks uses **23% less plastic** vs. standard flask design.



NEW:

Rounded shoulders provide:

- Improved scraping with fewer corners
- More efficient cell recovery
- Easier grip
- Better access to cap for removal and tightening
- Sterility assurance level (SAL) 10⁻⁶

Pipet Divot provides:

- No-tip liquid removal
- One-handed media addition and removal
- Access for 50mL pipets



These products included in the Corning EcoChoice™ program meet one or more of the following criteria:

- Recycled content (pre-consumer or post-consumer)
- Source reduction
- Intensification: products designed to enable greater cell production in a smaller footprint, resulting in less plastic per unit of output
- Manufactured using environmental attribute certificates (EACs)

To learn more about the EcoChoice program visit www.corning.com/EcoChoice.




Cat. No.	Surface (cm ²)	Neck	Cap	Qty/Cs
CON430641U-100ea	75*	Canted	Vented	100
CON430825-50ea	150	Canted	Vented	50
CON431080-50ea	175	Angled	Vented	50

*Pipet divot

Perfectly match with:

Corning® Cell Scrapers



- Useful for the manual harvesting of cells
- Blade design minimizes cell damage and ensures even contact with the growth surface
- Cell scrapers designed for use in flasks
- Individually wrapped, sterile, non pyrogenic

Cat. No.	Description	Blade Length (cm)	Handle Length (cm)	Qty/Cs	List Price	Promo Price
CON3010-100ea	Cell scraper, small	1.8	25	100	CHF 330.30	CHF 124.95
CON3011-100ea	Cell scraper, large	3.0	39	100	CHF 428.10	CHF 144.90

Offer valid until 31.08.24 - Prices are net, excluding VAT, and this promotion can't be combined with any other ongoing discounts or offer.